



CENTRIFUGE FOR WRV

Model N4000



Description

Refrigerated powerful cooling unit

- Automatic recognition of rotors
- Free rotor speed up to 4500 r / min or 3859 xg
- Intuitive graphic display 2 buttons. Memory 10 programs.
- Programming in t / min or acceleration xg, simultaneous display of two values
- Adjustable Acceleration: maximum or progressive.
- Braking 3 positions: maximum, progressive or simple inertia of the rotor.
- Automatic lid opening and audible signal activated at the end of cycle.
- Motor brushless, maintenance free, extremely quiet
- Equipped: a free 4-position rotor, 4 round nacelles, 4 racks and 4 test cells WRV 3 sampling areas

PTA SPAIN

C/ Gurutzegi N° 16 Pol. Belartza
20018 SAN SEBASTIÁN
(Guipúzcoa) - Spain

PTA FRANCE

14, rue Evariste Galois / ZAC les Ruieres
38320 EYBENS - France

PTA ASIA

50 East Coast Road #02-73 Roxy Square
428769 SINGAPORE - Singapore