



BAUER MC NETT Fractionators 4

Code. X10.503.400

For measure the weighted average fibre length of a pulp using the classification method.

Fractionators with 4 tanks with different sieves and integrated vacuum pump.

Standards : TAPPI T233 / SCAN M 6 / PAPTAC C5V



PTA SPAIN

C/ Gurutzegi N° 16 Pol. Belartza
20018 SAN SEBASTIÁN
(Guipúzcoa) - Spain

PTA FRANCE

14, rue Evariste Galois / ZAC les Ruieres
38320 EYBENS - France

PTA ASIA

50 East Coast Road #02-73 Roxy Square
428769 SINGAPORE - Singapore



PULP Testing Equipment

PAPER TESTING ASSOCIATION

SPECIFICATION

- Sturdy frame complete made of stainless steel on which are mounted 4 tanks with built-on motors.
- Integrated control box with fully automated test cycle
- Automatic start-stop with switch for adjusting the sequence time
- Quick clamping device for easy mounting the sieves
- Built-in vacuum pump for sucking out the water of the tanks after fractionating
- Holders for the filter papers with quick clamping

DELIVERY CONTENT

- Bauer McNett fractionator with 4 tanks and 4 sieves
- Sieves mounted on a sturdy sieve frame with 16, 30, 50 and 100 mesh
- Storage container for the sieves when not in use
- Water jet tubing hose for washing the containers during emptying
- Calibration certificate
- Operation manual
- Connection cable

OPTION :

- As option: sieve with 200 mesh

PTA SPAIN

C/ Gurutzegi N° 16 Pol. Belartza
20018 SAN SEBASTIÁN
(Guipúzcoa) - Spain

PTA FRANCE

14, rue Evariste Galois / ZAC les Ruires
38320 EYBENS - France

PTA ASIA

50 East Coast Road #02-73 Roxy Square
428769 SINGAPORE - Singapore

info@ptagroup.com

www.ptagroup.com



PAPER TESTING ASSOCIATION

PULP Testing Equipment



CONNECTIONS :

Electric : 230 V / 50 or also 110 V / 50 Hz or 60 Hz

Water : Drinking water quality min. 2 bar (200 kPa),
Drainage 50 mm diameter

WEIGHT AND DIMENSIONS (ONLY MACHINE):

1700 x 600 x 1600 mm (WxLxH) / 250 Kg

WEIGHT AND DIMENSIONS WITH PACKING BOX :

1900 x 800 x 1800 mm (WxLxH) / 350 Kg

PTA SPAIN

C/ Gurutzegi N° 16 Pol. Belartza
20018 SAN SEBASTIÁN
(Guipúzcoa) - Spain

PTA FRANCE

14, rue Evariste Galois / ZAC les Ruieres
38320 EYBENS - France

PTA ASIA

50 East Coast Road #02-73 Roxy Square
428769 SINGAPORE - Singapore

info@ptagroup.com

www.ptagroup.com