

Bauer McNett fibers classifier 5-unit | 070314

Compliant with standards

TAPPI T233, SCAN M6, PAPTAC C5V



Description

Apparatus for determining the fiber length distribution in pulps.

- Control screen with automated cycles.
- Easy and quick replacement of the grids.
- Stainless steel models with 3, 4 or 5 fractioning stations.
- Robust modular frame made of stainless steel with modular sorting units.
- Each sorting station is equipped with its drive motor.
- Built-in flowmeter.
- Vacuum pump included in the standard supply.
- Classification stations with quick opening system.
- Screening mesh for 3, 4 or 5 sorting stations (ASTM 16, 30, 50, 100 and 200 Mesh).

Characteristics

Contents of the delivery