


Bending stiffness tester | ST-10

Compliant with standards

ISO 2493,SCAN P29,DIN 53121,BS 3748,TAPPI T556,NFQ 03-048



Description

Device to determine the bending strength and the bending stiffness of paper and paperboard.

- Pneumatic tighten.
- Screening test length from 1-50 mm.
- Pre-settable bending angle (0.1 to 90 degrees).
- Angle adjustable speed from 0.1 to 10 ° / second.
- Load cell 10N- Range 10-10000mN- Resolution 0,5%.
- Large graphic touch screen.
- Clear acrylic security cover.

Characteristics

Electrical : 110V-230V, 50-60Hz

Contents of the delivery