


## Crush Tester | SC-500

### Compliant with standards

ISO 3035, ISO 3037, TAPPI, T808, T809, T811, T818, T821, T825, T839, FEFCO N° 6-8 & 11, AFNOR 23035, etc.



### Description

Device to carry out different type of compression tests on paper and corrugated board : CMT, RCT, CCT, ECT, FCT, PAT.

- Compression plates: 150 x 150 mm.
- Maximum force: 5000 N.
- Maximum separation between plates: 150 mm.
- User adjustable stroke limits.
- Safety stop button.
- Test speed between 1 and 200 mm/min.
- Automatic approach and return speed.

#### Characteristics

230 V / 50 Hz

#### Contents of the delivery

Crush Tester  
Electrical connection cable  
User manual