

Disintegrator | N6200_A

Compliant with standards

ISO 5263-1,SCAN C18 M2,TAPPI T205,PAPTAC C.6



Description

Apparatus for the disintegration of pulp in suspension.

- 3 safety systems.
- Digital revolution counter and digital display.
- Speed of revolution according to ISO and TAPPI.
- Top opening.
- Disintegrating pot in acrylic glass or stainless steel.

Characteristics

Electricity: 230 V / 50 Hz AC

Contents of the delivery

1 disintegrator
1 stainless steel or plastic pot