

## Formette Rapid Köthen | 0708A1



### Description

**Apparatus for the preparation of a laboratory standard hand sheet.**

- Touch screen to manage all the functions of the equipment.
- Automatic process controlled by PLC or manually.
- Polypropylene and stainless steel materials.
- Lightweight dryer cover for easy handling.
- Drying temperature of  $93 \pm 1^\circ\text{C}$  with white water recirculation according to standards
- Four models available: with 1, 2, 3 and 4 dryers.

## Formette Rapid Köthen | 0708A1

### Technical specifications

The Rapid Köthen is used for the production of standardised paper sheets for the laboratory.

The paper sheets are used to determine the physical and optical properties of the pulp.

The formette is equipped with a 205 mm diameter formation column with dryers. The production of the handsheet does not require any operator intervention.

**Two models are available: manual and automatic.**

#### Characteristics

#### Contents of the delivery