


Laboratory pulper V.A 8 liters | N2101-A08



Description

Device for the defibration of dry samples in laboratory scale.

- Different volumes available: 5/8/17/45 liters corresponding to useful volumes of 3/5/10/30 liters.
- Made of stainless steel.
- Fast defibration of virgin and secondary fibres.
- For a concentration between 2 and 7% dry matter.
- Variable turbine speed from 50 to 1500 rpm.
- Programmable automatic shut-off timer.
- Acrylic lid to avoid splashing.

Characteristics

110 – 230 V / 50 – 60 Hz

Contents of the delivery