

## SCT apparatus | SCT-07

### Compliant with standards

ISO 9895, TAPPI T826, DIN 54518, SCAN P46, BS 7325



### Description

Device to measure the short and zero-span compression (SCT) on paper and board stripes.

- Ratio calculation between MD and CD samples.
- Range: 1 \_ 500 n equal to 0,06 \_ 33,3 kN/m.
- Sample width: 15 mm.
- Span of clamps: 0,7 mm.
- Clamping pressure: 2300 +/- 50 N (adjustable).

#### Characteristics

110 – 230 V / 50 – 60 Hz  
600 kPa

#### Contents of the delivery

Short Span Compression Tester  
Connection cable  
User manual