

Somerville shive content analyser | N7100-A

Compliant with standards

TAPPI T275, TAPPI UM242, PAPTAC C11.P



Description

Apparatus for the determination of shive content in pulp.

- Stainless steel and POM construction.
- Different sizes of grids are available.
- Control (on/off/water), timer for 20 min sieving time (motor and water).
- Thermal overload protection integrated in the electrical control unit (IP65).
- Air gun supplied complete with tube and self-sealing connection for cleaning.

Characteristics

Contents of the delivery

Device delivered with feet
Calibrated orifice
Electric box
Hand shower and hoses