

## Valley type beater | 070113

### Compliant with standards

ISO 5264/1, SCAN C25, TAPPI T200m, T205m,



### Description

Traditional Dutch pile type laboratory refiner according to the standard.

- Made entirely of stainless steel.
- Refiner roller speed: 500 r.p.m.
- Tub volume 23 liters.
- Maximum sample capacity 360 g.
  
- Refining tank and structure made entirely of stainless steel.
- Good reproducibility of refining pressure, adjustable by weights.
- Tub volume 35 liters. Refiner roller speed: 500 r.p.m.
- It allows to calculate the power used to obtain the desired refining degree.

#### Characteristics

#### Contents of the delivery