

Laboratory Disc Refiner PTA | 070111



Description

A disc machine for refining pulp similar to the refiners used in production.

- Ø Disc : 200mm - 8inches.
- Speed: 3.600 r.p.m.
- Stainless steel tank of 45L (other volumes on request).
- Disc refiner for pulp in MP process conditions.
- Consistency of the pulp between: 2-8%.
- Ø disc: 200mm - 8inches.
- Speed: 3.600 r.p.m.
- 45L Stainless steel tank (other volumes on request).
- Stainless steel disc chamber and holders.

Characteristics

Contents of the delivery